

F&P Airvo 3

Optiflow high flow therapy device

F&P Airvo™ 3

Optiflow™ High Flow System


Next generation humidifier
with integrated flow generator

www.fphcare.com/airvo3

Fisher & Paykel
HEALTHCARE

F&P Airvo 3 Optiflow High Flow System.

A humidifier with integrated flow generator, delivering respiratory support and airway hydration with a focus on mobility and continuity across the patient journey.



Standardised consumables.

Compatible with existing Airvo 2 patient consumables, facilitating integration into an existing Airvo 2 fleet.

Extended compatibility to smaller Optiflow Junior 2 nasal interfaces

Simple touch screen display.

Quick selection and adjustment of therapy settings.

Visualise patient trends and graphs.



Single Optiflow mode.

One mode to enable Optiflow high flow applications across adults and pediatric patients.

Flow range extended from 2 to 70 L/min to treat a wider group of patients.*

Flexible oxygen inputs.

High-pressure oxygen input enables precise FiO_2 adjustment directly on the device, automatically adjusting oxygen input to meet target FiO_2 .

Alternative low-pressure oxygen input also available.



Integrated battery.

Enables 40+ minutes of operating time to facilitate transport and patient mobilisation throughout the hospital environment.



Download
the Airvo 3
Simulator App



Hardware and accessories	Product Code
Airvo 3 device	PT301XX*
Mobile pole stand	900PT421
Storage basket	900PT42E
Oxygen-bottle holder regular/large	900PT423/L
HPO Dual-input manifold	900PT460D

* Indicates country specific model. ^ Indicates a specific size of interface.

Patient consumables	Product Code
AirSpiral™ tube and chamber kit	900PT561
AirSpiral tube and chamber kit including Airvo neb adapter	900PT562
Optiflow+ nasal cannula	OPT94X*
Optiflow+ Tracheostomy direct connect	OPT970
Optiflow Junior 2 nasal cannula (M, L, XL)	OJR414/416/418

F&P, Airvo, AirSpiral, Duet and Optiflow are trademarks of Fisher & Paykel Healthcare Ltd.
For patent information, see www.fphcare.com/ip.
629872 REV A © 2023 Fisher & Paykel Healthcare Limited. This document is not intended for USA.